

Validation of the Analytical Sensitivity of Commercial Scale Lots of the Verax PGD*prime*® Test for Bacteria in Platelets

Lisa Shinefeld, David LaVerda, Nancy Best, Johnny Lisitu,
Gary Tambolleo, Remo Vallejo

Verax Biomedical Incorporated, Marlborough, MA


Validation of the Analytical Sensitivity of Commercial Scale Lots of the Verax PGDprime® Test for Bacteria in Platelets

Lisa Shinefeld, David LaVerda, Nancy Best, Johny Lisitu, Gary Tambolleo, Remo Vallejo

Verax Biomedical Incorporated, Marlborough, MA

Background

- The PGDprime rapid test for bacteria in platelets was developed as an improvement of the PGD® test currently used as a Safety Measure for platelet transfusion.
- It was designed to provide a simpler test procedure compared to the current product and improved specificity through the use of intact IgG for capture and F(ab')₂ fragments for detector antibodies for 5 out of 6 test lines.
- Using 3 sequential cGMP lots, the Analytical Sensitivity of the PGDprime test was compared to that of PGD.



PGDprime Test

- PGD Procedure

 - Add Reagent 1 (lysing agent) to platelet sample
 - Centrifuge
 - Decant supernatant
 - Add Reagent 2 (base) to pellet
 - Disrupt pellet and mix
 - Add Reagent 3 (neutralizer) and vortex
 - Transfer to PGD test device.
- PGDprime Procedure

 - Add Reagent 1A (base) to platelet sample. Invert to mix
 - Add Reagent 1B (neutralizer) to sample. Invert to mix.
 - Transfer to test device.
 - Add Chase buffer (Reagent 2).

No centrifuge is required.

Methods

Phase 1: Sensitivity in Apheresis Platelets

- A blinded, frozen 32-member bacterial test panel was prepared. The panel comprised 3 levels of each PGD bacterial claim strain so that, after reconstitution in platelet samples, one member was below the limit of detection (LoD) of the PGD assay (“Lower”), one member was at or slightly above the LoD of the assay (“Middle”), and one was above the LoD of the PGD assay (“Upper”). Two negative samples were included.
- Panel members were thawed and added 1:21 to 10 in-date apheresis platelets in plasma.
- The reconstituted panel members were tested on one lot of PGD and on 3 lots of PGDprime. Bacteria included *B. cereus*, *C. perfringens*, *S. aureus*, *S. epidermidis*, *S. agalactiae*, *E. aerogenes*, *E. coli*, *K. pneumoniae*, *P. aeruginosa*, and *S. marcescens*.

Bacterial Test Panels

Bacterial ID	Level	PGD LoD Claim CFU/mL	Panel CFU/mL	Log ratio to LoD
<i>Bacillus cereus</i>	Upper	1.2E+04	7.8E+04	0.8
	Middle		2.7E+04	0.4
	Lower		3.4E+03	-0.5
<i>Clostridium perfringens</i>	Upper	8.9E+04	7.8E+05	0.9
	Middle		2.4E+05	0.4
	Lower		3.6E+04	-0.4
<i>Staphylococcus aureus</i>	Upper	8.2E+03	7.7E+04	1.0
	Middle		1.8E+04	0.3
	Lower		2.1E+03	-0.6
<i>Staphylococcus epidermidis</i>	Upper	9.2E+03	4.9E+04	0.7
	Middle		2.7E+04	0.5
	Lower		1.9E+03	-0.7
<i>Streptococcus agalactiae</i>	Upper	5.5E+04	5.1E+05	1.0
	Middle		1.6E+05	0.5
	Lower		1.8E+04	-0.5
<i>Enterobacter aerogenes</i>	Upper	1.0E+04	5.6E+04	0.8
	Middle		3.3E+04	0.5
	Lower		1.1E+03	-0.9
<i>Escherichia coli</i>	Upper	2.8E+04	1.2E+05	0.6
	Middle		5.6E+04	0.3
	Lower		4.2E+03	-0.8
<i>Klebsiella pneumoniae</i>	Upper	2.0E+04	1.8E+05	1.0
	Middle		6.1E+04	0.5
	Lower		2.4E+03	-0.9
<i>Pseudomonas aeruginosa</i>	Upper	8.2E+03	3.0E+04	0.6
	Middle		1.7E+04	0.3
	Lower		2.6E+03	-0.5
<i>Serratia marcescens</i>	Upper	8.6E+05	7.6E+06	0.9
	Middle		2.5E+06	0.5
	Lower		1.7E+05	-0.7
Negative Control	Neg	0.0E+00	0.0E+00	0
	Neg		0.0E+00	0

Phase 2: Sensitivity in Other Platelet Types

An LoD Challenge Study in other platelet types was conducted by spiking the mid-level panel members into 6 units each of the following platelet types:

- pre-storage pools of leukocyte-reduced random donor platelets (PSP)
- post-storage pools of non-LR random donor platelets (NLR pools)
- apheresis platelets in platelet additive solution (PAS)
- An apheresis platelet unit was also spiked with the mid-level panel members as a control.
- S. oralis* was included as a new detection target for PGDprime. It was tested at a level typically associated with previously reported serious septic events with platelets.

Results/Findings

Phase 1: Apheresis Platelets Spiked with Bacterial Panels

LoD Comparison (Reactive Results/Units Tested)			
Bacteria	Level (mid=LoD)	PGD	PGDprime (3 Lots)
<i>B. cereus</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	26 / 30
<i>S. aureus</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	22 / 30
<i>S. epidermidis</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	30 / 30
<i>S. agalactiae</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	30 / 30
<i>C. perfringens</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	1 / 10	29 / 30
<i>E. aerogenes</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	2 / 30
<i>P. aeruginosa</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	4 / 30
<i>E. coli</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	2 / 30
<i>S. marcescens</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	1 / 10	30 / 30
<i>K. pneumoniae</i>	Upper	10 / 10	30 / 30
	Middle	10 / 10	30 / 30
	Lower	0 / 10	0 / 30
Negative		0 / 20	0 / 60*

*2 initial reactivities, no repeat reactivities

Phase 2: LoD Challenge in Other Platelet Types

Bacteria	Platelet Type Tested	PGDprime (3 Lots)
<i>B. cereus</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>S. aureus</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>S. epidermidis</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>S. agalactiae</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>C. perfringens</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>E. aerogenes</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>P. aeruginosa</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>E. coli</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>S. marcescens</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>K. pneumoniae</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
<i>S. oralis</i>	PSP	30/30
	NLR pool	18/18
	apheresis in PAS	21/21
Negative	PSP	0/20
	NLR pool	0/18
	apheresis in PAS	0/21

Conclusions

- All samples at and above the LoD of PGD were detected by PGD and by all 3 lots of PGDprime in apheresis platelets.
- For some strains, PGDprime detected bacterial concentrations lower than the LoD of PGD.
- In the LoD challenge study, PGDprime detected all species at the LoD of PGD in all platelet types.

